

XINEOS

Dynamic CMOS X-Ray Flat Detectors
for Medical Radiology



Xineos - CMOS Flat X-Ray Detectors

IMPROVED MEDICAL DIAGNOSTIC AT A FRACTION OF THE DOSE

Teledyne DALSA Xineos offers a complete portfolio of CMOS flat X-Ray detectors tailored specifically to meet the demanding needs of dynamic diagnostic and interventional medical applications service, support and manufacturing.

- General Surgery
- Cardiology
- Vascular procedures
- Gastro-intestinal Fluoroscopy
- Traumatology
- Orthopedics
- Electrophysiology
- Urology



CONVENIENT SQUARE FIELD-OF-VIEW – NEXT GENERATION X-RAY TECHNOLOGY

The Xineos flat detector family spans the popular square field-of-view dimensions.

- 13x13cm²
- 15x15cm²
- 22x22cm²
- 30x30cm²

Replaces the incumbent technology Image Intensified CCD (IICCD) systems with circular diameter of 6", 9" or 12".

Evolution of the Dynamic X-ray Detection

Image Intensified CCD
++ Image quality at Low dose
– Bulky form factor
– Distortions



1948

TFT Flat detectors
– Image quality at Low dose
++ Slim form factor
++ No Distortions




1997

TeledyneDalsa's CMOS Flat detectors
++ Image quality at Low dose
++ Slim form factor
++ No Distortions



2010

Xineos – Family of Interventional CMOS Detectors




EASY TO DEPLOY – VERSATILE FUNCTIONALITY

- No protruding connectors due to connector bay
- Slim and elegant design
- EMC resistant
- Switchable saturation dose
- Seamless and instant readout mode switching
- Sensitivity independent of pixel binning mode
- Low weight
- Flexible zoom modes
- CE and UL certifications
- Software libraries for interfacing and real-time image pre-processing
- Stable offset and gain calibration, no temperature dependency



Xineos – True New Digital X-ray Advancements

The Xineos family of CMOS X-Ray detectors is the first commercially available detectors that combine the low dose image quality (IQ) of the IICCDs and the multiple advantages of amorphous-silicon (aSi) flat detectors, including a slim form factor, no geometric distortions, and a larger field-of-view (FOV).



- **LOW DOSE AND HIGH IQ AT WIDE DOSE SPECTRUM**

The Xineos Family of medical detectors features the highest Detective Quantum Efficiency (DQE) at very low dose levels, surpassing the previous ‘golden standard’ image quality performance of IICCDs. The very low noise of the CMOS material and the proprietary active pixel architecture of Teledyne DALSA’s CMOS detectors assure superior signal-to-noise ratio (SNR) at low dose levels compared to a-Si-based and even other CMOS-based competitors.

- **HIGH DYNAMIC RANGE**

Utilizing the latest generation proprietary radiation-hard active pixel design with selectable pixel saturation charge capacitors, the Xineos detectors are extremely versatile and suitable for all types of interventional applications using just a single detector. Those detectors can seamlessly change from high sensitivity to high dynamic range modes, without compromising on image quality or patient dose, and without elaborate and time consuming calibration routines.

- **REAL-TIME IMAGING AT FULL RESOLUTION AND FULL FIELD-OF-VIEW**

Enabled by high speed CMOS integrated circuits, the Xineos family of dynamic X-ray detectors sets the industry benchmark for speed at full field-of-view (FOV). All detectors of the family support 30fps at full resolution and full FOV. The images of the Xineos detectors are always lag- and artifact-free, due to the high electron mobility of the crystalline silicon material, which is not the case when using any other X-Ray imaging technology. Detector’s offset and gain calibration is very stable over the operating temperature range.

Improved Return on Investments

- **HIGH RESOLUTION**

High spatial resolution is achieved in the Xineos family of detectors, offering a small pixel pitch of 99µm across all FOV models. Teledyne DALSA's advanced active pixel design is responsible for the very high 85% pixel fill factor, even at such small pixel sizes. The small pixel pitch together with a proprietary scintillator assembly process also contributes to high MTF performance of the detectors.

- **LOW DISSIPATION POWER**

The high level of component integration and the low resistivity of the crystalline CMOS material in the Xineos family detectors contribute to the very low power dissipation, which makes any kind of active cooling in the medical systems obsolete. This saves high running and maintenance costs for expensive liquid cooling systems, as well as improves up-time and reliability of critical imaging equipment.

- **INSTANT START-UP TIME AND STABLE OPERATION**

The extremely low dark current of the Xineos detector family is responsible for a very stable operation of the detectors right from the start, where the effects of system warming up does not influence the operation of the detector. Teledyne DALSA's detectors are fully operationally running 30 seconds after switching on the power, typically faster than the system host boot time.

- **LONG LIFETIME AND COST REDUCTION**

The Xineos detector family utilizes a high integration level of discrete peripheral electronic components into the pixel array. This drastically reduces the number of components and interconnects, which in turn reduces the detector cost and the assembly complexity while at the same time improving the product reliability. The proprietary active pixel of the Xineos detector family is radiation-hard by design, meaning that it is resistant to performance degradation caused by X-Ray radiation. This enables long operating lifetime and less frequent calibration routines.



Improved Patient Care, Optimized Diagnostics Simplified Deployment



RADIOLOGISTS

- High image quality at very low dose – lower dose for patient and staff, higher success rate of procedures, higher patient throughput
- System operational at under 30 seconds – immediate start of the system – higher success rate of emergency procedures
- Reliable and stable detector operation – no calibration routines and loss of time, higher throughput
- Seamless mode switching – no loss of time, easy system operation



HOSPITAL ADMINISTRATORS

- High image quality at very low dose –
 - » physicians can run more procedures for the same dose quote per annum
 - » better hospital positioning due to higher success rate of procedures
 - » better ROI thanks to higher patient throughput
- Reliable and stable detector operation – no calibration, satisfied physicians, faster procedures
- Seamless mode switching – more and faster procedures
- Total cost of ownership minimization – longer detector lifetime, less service costs for the system



ORIGINAL EQUIPMENT MANUFACTURERS

- Added value to the system due to:
 - » Excellent IQ at low dose performance
 - » Superior spatial resolution
 - » Faster start up time
 - » Seamless mode switching
- Faster and easier detector integration
 - » Real-time image preprocessing, without external processing box
 - » No active cooling required
- Reliable systems, less maintenance, fewer support calls
- Possibilities for new applications due to high speed capabilities
- Local design-in and aftersales support

XINEOS

FAMILY DATA SPECIFICATIONS

PARAMETER	UNIT	XINEOS-1313FL	XINEOS-1515	XINEOS-2022HR	XINEOS-2222HS	XINEOS-3030HR	XINEOS-3030HS
GENERAL							
TECHNOLOGY		CMOS Active Pixel	CMOS Active Pixel	CMOS Active Pixel	CMOS Active Pixel	CMOS Active Pixel	CMOS Active Pixel
PIXEL PITCH	[um]	100.1	99	99	151.8	99	151.8
PIXEL CAPACITY MODES	[#]	1	2	2	3	2	3
ACTIVE AREA	[mm]	131x131	153x153	206x224	215x215	296x296	296x296
RESOLUTION	[px]	1316x1312	1548x1548	2065x2237	1416x1420	2994x2997	1952x1952
IMAGE QUALITY							
SATURATION DOSE, RQA5 (per mode)	[uGy]	2	3 / 14	5 / 23	2 / 10 / 45	5 / 23	2 / 10 / 45
MTF @ 1LP/MM	[%]	60%	60%	60%	60%	58%	58%
DQE @ 0LP/MM, RQA5	[%]	70%	70%	70%	70%	66%	66%
RANDOM NOISE (HFW)	[DN]	6	5 / 3	5 / 3	5 / 3 / 2.5	5 / 3	5 / 3 / 2.5
DYNAMIC RANGE (HFW)	[dB]	69	71 / 74	71 / 75	70 / 75 / 77	71 / 75	70 / 75 / 77
IMAGE LAG, FIRST FRAME @ 30fps	[%]	0.1%	0.1%	0.1%	0.1%	0.1%	0.1%
BANDWIDTH							
DATA INTERFACE	[-]	GigE / CameraLink	GigE / CameraLink	CameraLink	GigE	CameraLink	GigE
IMAGE BIT DEPTH	[bits]	14	14	14	16	14	16
FRAME RATE - 1x1 (CamLink)	[fps]	45	30	31	–	31	–
FRAME RATE - 2x2 (CamLink)	[fps]	60	60	56	–	63	–
FRAME RATE - 1x1 (GigE)	[fps]	30	22	–	30	–	15
FRAME RATE - 2x2 (GigE)	[fps]	60	60	–	90	–	60
ROI READOUT		(X, Y) flexible	(X, Y) flexible	(X, Y) flexible	(X, Y) flexible	(X, Y) flexible	(X, Y) flexible
POWER CONSUMPTION							
POWER SUPPLY	[Vdc]	12	11..26	11..26	11..26	11..26	11..26
POWER CONSUMPTION	[W]	11	8	15	15	18	18
ACTIVE COOLING	[-]	NO	NO	NO	NO	NO	NO
INTEGRATION							
FOOTPRINT (WxH)	[mm]	188x150	224x176	292x235	292x235	377x327	377x327
THICKNESS (IN/OUTSIDE IO)	[mm]	58 / 35	45 / 36	56	56	58	58
WEIGHT	[kg]	2.5	3.0	6	6	9	9
EXTERNAL INTERFACE MODULE	[Y/N]	NO	NO	NO	NO	NO	NO
ENVIRONMENTAL							
CERTIFICATION		–	CE/UL	CE/UL	CE/UL	CE/UL	CE/UL
EMC (IEC 60601-1-2)		CLASS B	CLASS B	CLASS B	CLASS B	CLASS B	CLASS B
OPERATING TEMPERATURE	[°C]	+10..+40	+10..+40	+10..+40	+10..+40	+10..+40	+10..+40
STORAGE TEMPERATURE	[°C]	-10..+55	-10..+55	-10..+55	-10..+55	-10..+55	-10..+55
HUMIDITY (NON-CONDENSING)	[% R.H.]	20..80	20..80	20..80	20..80	20..80	20..80
X-RAY ENERGY RANGE	[kV]	40..120	40..125	15..150	15..150	15..150	15..150

World-Class Capability...

Design, Engineering and Manufacture

LEADING PROVIDER AND PREFERRED PARTNER

Image Sensor Product Solutions in Professional- and Mid-Range segments

- Medical/Dental X-Ray
- Non-Destructive Testing
- Industrial & Scientific X-Ray
- Human Vision

IMAGING HERITAGE OF 35+ YEARS

Design, Development & Manufacturing of

- X-Ray Detectors
- CCD & CMOS Image Sensors
- Chipset & Application Reference Designs
- Customer Application Support

A QUALITY CERTIFIED COMPANY

- ISO 9001-2008
- ISO 13485
- ISO 14971
- IEC 60601-1 3rd edition
- UL & CE certification
- RoHS compliancy

● **Headquarters + Manufacturing** ● **Manufacturing** ● **Sales & Customer Support**



www.teledynedalsa.com

For more detailed product information please visit:

www.teledynedalsa.com/ls or contact us at sales.sensors@teledynedalsa.com

Teledyne DALSA has its corporate offices in Waterloo, Canada

Teledyne DALSA reserves the right to make changes at any time without notice. Teledyne DALSA © 2017.

AMERICAS

3400 Garrett Drive
Santa Clara, CA
USA 95054

Tel: 408-736-6000
Fax: 408-736-6001

EUROPE

High Tech Campus 27
5656 AE Eindhoven
The Netherlands

Tel: +31 40 2599000
Fax: +31 40 2599005

ASIA PACIFIC

Ikebukuro East 13F
3-4-3 Higashi Ikebukuro
Toshima-ku, Tokyo
Japan 170-0013

Tel: +81 3-5960-6353
Fax: +81 3-5960-6354

Shanghai Industrial Investment Building
Room G, 20F, 18 North Cao Xi Road
Shanghai, China 200030

Tel: +86-21-64279081
Fax: +86-21-64699430